

Gas Analysis | HPR-20 R&D  
Application Note AN-10011.1


## NO<sub>x</sub> Detection

The detection of NO and NO<sub>2</sub> using the HPR-20

### Summary

The Hiden HPR-20 gas analysis system has been used in a variety of applications in the fields of catalysis and thermal analysis. One important aspect of catalysis research is the reduction of NO<sub>x</sub> in such applications as car exhausts and power plants. This research involves many challenges in meeting current emission legislation whilst keeping catalyst cost low. Several NO<sub>x</sub> reduction technologies have emerged in previous years, some of which have now been commercialized including selective catalytic reduction (SCR) with NH<sub>3</sub>, urea and hydrocarbons, and NO<sub>x</sub> storage and reduction (NSR). NSR technology consists of two cyclic steps that occur on a lean NO<sub>x</sub> trap (LNT) catalyst. The LNT catalyst readily stores NO<sub>2</sub> as compared to NO [1,2]. For this reason, NO should be oxidized to NO<sub>2</sub> to achieve an acceptable level of NO<sub>x</sub> storage.

This application note describes the measurement of NO<sub>2</sub> using the m/z 46 peak and the NO/NO<sub>2</sub> ratio in a research application using the Hiden HPR-20. The data highlights the ability of the QIC series gas analysis systems to deconvolute and quantify low levels of NO in high concentrations of NO<sub>2</sub>.



**Figure 1: Profile scan for 10000ppm NO<sub>2</sub>.**

## Introduction

The Hiden HPR-20 gas analysis system is configured for continuous analysis of gases and vapours.

The Hiden QIC quartz-lined sampling interface operating up to 200°C provides fast response times of less than 300 millisecond for most common gases and vapours, including water vapour. For this particular application note, the HPR-20 system had the following specifications:


- i) QIC capillary inlet with 2 metre heated capillary, operating up to 160°C
- ii) Quadrupole probe with 100 amu mass range capability.
- iii) Triple-stage mass filter
- iv) Corrosion resistant pumping system

## Test Data

Figure 1 shows the result of measuring NO<sub>2</sub> in a carrier gas of He. It clearly demonstrates that a distinct peak at m/z 46 can be measured using the HPR-20 system.


This is important as it is known that NO<sub>2</sub> rapidly and easily fragments to NO<sup>+</sup> making analysis of any NO<sub>x</sub> mixtures containing NO and NO<sub>2</sub> difficult to deconvolute thereby increasing the uncertainty in the NO and NO<sub>2</sub> concentration results. Assuming the major peak at m/z 30 represents 100%, the peak at m/z 46 is equivalent to 15% giving a 46/30 ratio of 1:6.6.

Another potential problem of the facile fragmentation of NO<sub>2</sub> to NO<sup>+</sup> is the measurement of low concentrations of NO in a high background of NO<sub>2</sub> (a vital measurement in developing NO<sub>x</sub> reduction technologies such as the LNT catalyst). The NO<sub>2</sub> produces a high background at m/z 30, which could potentially mask any signal due to the NO in the gas mixture. Figure 2 shows an example MID scan for a range of different NO concentrations measured at m/z 30 with a background level of 1000ppm NO<sub>2</sub>. It can clearly be seen that, for the Hiden HPR-20 system, detection limits better than 12ppm are achievable.



**Figure 2: Detection Limit of NO in 1000ppm NO<sub>2</sub> (Inset: Expansion of detection of 12ppm NO in 1000ppm NO<sub>2</sub>). Primary Axis represents Partial Pressure at m/z 30, Secondary Axis represents Partial Pressure at m/z 46.**

Correlating the NO peak signal intensity with the NO concentration level provides quantitative information on the NO detection levels as shown in Figure 3. The high degree of linearity indicates a high level of confidence in the detection of ppm levels of NO in a high background of NO<sub>2</sub>.



**Figure 3: Relationship between partial pressure for m/z 30 and NO concentration (Note: The background at m/z 30 due to NO<sub>2</sub> has been subtracted).**

## Conclusions

- The Hiden HPR-20 shows high sensitivities for NO<sub>2</sub> detection using the m/z 46 peak despite a large degree of fragmentation.
- Low levels of NO detection in a high background of NO<sub>2</sub> can be achieved without the requirement for other additional or tandem techniques.
- The specific data as detailed in this application note shows detection levels better than 12ppm of NO are attainable – a significant benefit of the Hiden QMS system.

## References

- [1] S. Erkfeldt, E. Jobson, M. Larrson, *Top. Catal.*, 2001, 16/17, 127.
- [2] W.S. Epling, L.E. Campbell, A. Yezerets, N.W. Currier, J.E. Parks II, *Catal. Rev.*, 2004, 46, 164